

DIA-PROT



- Designed to be used to complement Protaper Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#F1	ML 150-S601
#F2	ML 150-S602
#F3	ML 150-S603
#F4	ML 150-S604
#F5	ML 150-S605
#F1/F3	ML 150-S691
#F4/F5	ML 150-S692
#F1/F5	ML 150-S693

(Pkg. of 60)

DIA-PRO R



- Designed to be used to complement RECIPROC® Files.
- Millimeter Marked Points

Size	ITEM No.
#25	ML164-605
#40	ML164-608
#50	ML164-610
#Asst(25/40/50)	ML164-691

(Pkg. of 60)

DIA-PRO W



- Designed to be used to complement Wave One Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	ML165-601
#Primary	ML165-602
#Large	ML165-603
#Asst(Small/Primary/Large)	ML165-691

(Pkg. of 60)

DIA-PINK



- Shorter than usual to suit slightly shorter canals
- available in 20 mm length
- Widely used in Europe
- One hand slide open lid

Size	ITEM No.
#XXF	103-901
#XF	103-902
#F	103-903
#M	103-904
#L	103-905
#XXF/M	103-991

(Pkg. of 100)

FEATHERED TIP



- Designed to be used as auxiliary points to the master cone technique
- Extra fine taper
- Hand Jig Rolled

Size	ITEM No.
#XF	126-601
#FF	126-602
#MF	126-603
#F	126-604
#FM	126-605
#M	126-606
#ML	126-607
#L	126-608
#XL	126-609
#XF/M	126-691

(Pkg. of 100)

DIA-GTP




- Designed to be used to complement Profile GT Greater Taper Files.
- Available in .06, .08, .10, .12, taper.

Size	ITEM No.
DIA-GTP .06 #20	112-601
DIA-GTP .08 #20	112-602
DIA-GTP .10 #20	112-603
DIA-GTP .12 #20	112-604
DIA-GTP .06/.12 #20	112-691

(Pkg. of 60)

DIA-GP GAUGE

- **ITEM #781-001** 
- Designed to size a perfect master cone from any Gutta Percha Points.
- Drop it into the apical size and cut off with a blade.
- Autoclavable
- Dimensions: 92 x 3 x 20m/m(W x H x D)

