

DiaCem™ TempNE

NON EUGENOL TEMPORARY CEMENT

NEW

- Non-Eugenol Temporary cement for trial cementing restorations or temporary crowns and bridges



Features

- Easy and precise direct dosing
- No hand mixing and Stable mix of base and catalyst viscosity with the automix tip
- High bond strength to prevent leakage and sensitivity to cold and heat
- Easy removal for the patient comfort
- Low film thickness with high compressive and flexural strength
- Fast working time(3 min) and setting time(5 min in the mouth)

Packages

- ITEM#2004-2101 2 Syringes of 10g, 10 Automix tips
- ITEM#2001-6112 20 Automix tips

DIAPASTE™

CALCIUM HYDROXIDE PASTE WITH BARIUM SULFATE



Characteristics

- Premixed Calcium Hydroxide + Barium Sulfate paste for root canal treatment.
- Temporary root filling.
- Root canal filling for primary teeth.
- Use after pulp extirpation treatment.
- Filling on the pulpotomy.
- Temporary Pulp Capping.

Benefits

- Water Soluble: Easy to clean and remove, reduce operative time
- Easy and quick access
- Radiopaque
- Antibacterial

Packages

- ITEM #1001- 402 Regular Kit
- 2g of paste and 10 disposable tips with cap
- ITEM #1001- 315 Blue Disposable Syringe Tips
- 20 disposable tips in a blister package

Tip Specification

