

Endo-Clear™**WATER BASED EDTA ROOT CANAL CLEANSER**

Endo-Clear™ is the EDTA root canal cleanser, colorless and clear without foreign matters and effectively cleans the canal and removes smear layer. Chelating effect of Endo-clear™ softens root canal wall, dissolves calcium and removes it.

Characteristics

- Debridement in the root canal
- Dissolution of pulp remnants
- Removal of smear layer
- Microbicidal treatment within the canal
- Water based type

Benefits

- Chelating effect of Endo-Clear™ make root canal wall softened and calcium dissolved.
- Easy to enlarge the canal and clean it.
- Endo-Clear™ perfectly removes smear layer.
- Superior microbicidal treatment within the canal
- Excellent dissolution effect of pulp remnants
- Being based on water, more easier to get to apical area where hard for cream type to reach. So it effectively removes smear layer around apical area better than cream type.

Packages

- | | | |
|------------------------------------|------------------------------|------------------------------------|
| Regular Package | | |
| • ITEM # 2002-5001 Regular Package | EDTA 17% Solution 100 ml x 1 | Dia-Flex (=Irrigation syringe) x 1 |
| | Disposable tips x 3 | |
| • ITEM # 2002-5002 Refill Package | EDTA 17% Solution 100 ml x 1 | |

**DIAPREP™ Plus****EDTA ROOT CANAL PREPARATION CREAM**

EDTA Cream with Urea Peroxide, used as an irrigation agent for chemo-mechanical preparation of the root canal.

Characteristics

- Debridement in the root canal
- Dissolution of pulp remnants
- Removal of smear layer
- Facile movements of canal instruments
- Microbicidal treatment within the canal
- Clearance of dentin tubules and apical area

Benefits

- EDTA removes smear layer, softens dentin, and facilitates the removal of calcific obstruction.
- Urea Peroxide reacts with Sodium Hypochlorite to produce oxygen bubbles that flush out debris within the canal.
- The lubricant medium facilitates instrument movement within the canal to eliminate instrument binding or breakage.
- Unclogged dentin tubules and apical area permit better penetration of root canal medicaments and fillings.
- Provides rapid cleansing and enlarging of canals due to EDTA, which softens walls and smooth calcium deposits.
- Easier placement into canals than liquids.

Packages

- ITEM # 2002-1100 - 2 syringes of 6g, 10 disposable tips and 1 protective cap
- ITEM # 2001-6111 - 20 Disposable tips in a polybag

