

DIAFIL CORE™ AUTOMIX DUAL CURED CORE BUILD-UP COMPOSITE **NEW**

Dual-curing(Light and self-curing) core build-up and radiopaque two-component composite materials in Automix delivery system



Benefits

- Comfortable dispenser: Automix syringes deliver a correct and consistent mix on every uses
- Easy and precise application with Automix and Intraoral tips
- Dual curing: Light cure for 40 seconds or Chemical cure within 4 minutes
- High compressive strength composite resin and long-term durability
- Cuts and wears like dentin
- Superior bond strength with DiaPlus
- Available 2 colors : Natural color(A3 shade) for the optimal esthetic results and Blue color

Packages

- **ITEM#2001-8101 Natural Color(A3 shade)**
- **ITEM#2001-8102 Blue Color**
 - 2 syringes of 9g
 - 10 Intraoral tips
 - 10 Automix tips
- **ITEM #2001-8109**
 - 20 Automix tips and 20 Intraoral tips

DIAFIL CORE™ SELF-CURING CORE BUILD-UP COMPOSITE

Self-curing and long-term fluoride releasing composite materials consist of two-components, base plus catalyst for core build-up system



Benefits

- Provides good handling properties for easier use
- Offers more faster self cure set time
- Smooth and creamy materials help fill it into needle tips very easily
- High level of radiopacity

Packages

- **ITEM #2001-7101**
 - 10g Base + 10g Catalyst
 - 15 Spatulas
 - 50 Mixing pads
 - 10 Needle tips