

ENDODONTIC NICKEL TITANIUM HAND FILES/PLASTIC PACKAGE

Increased Flexibility

: NiTi Shaft is made of premium Nickel-Titanium with 5 times the flexibility compared to Stainless Steel Files.

Strength

: NiTi has 10 times the stress resistance compared to Stainless Steel Files.

Convenience

: Produces smooth and progressive enlargement of root canal.
NiTi follow existing canal path and has memory so there is no need to Pre-Bend the files.

Efficiency

: Color Coded handles and pre-inserted silicone stops.

Safety

: NiTi resists spontaneous breakage.
Increased angular rotation reduces fracture.
NiTi decreases zipping and stripping.



NITI REAMERS 
NICKEL-TITANIUM REAMERS(.02 TAPER)

size	21MM	25MM
#8	NT501-101	NT501-201
#10	NT501-102	NT501-202
#15	NT501-103	NT501-203
#20	NT501-104	NT501-204
#25	NT501-105	NT501-205
#30	NT501-106	NT501-206
#35	NT501-107	NT501-207
#40	NT501-108	NT501-208
#45	NT501-109	NT501-209
#50	NT501-110	NT501-210
#55	NT501-111	NT501-211
#60	NT501-112	NT501-212
#70	NT501-113	NT501-213
#80	NT501-114	NT501-214
#15/40	NT501-191	NT501-291
#45/80	NT501-192	NT501-292

(Pkg. of 6)

NITI K-FILES 
NICKEL-TITANIUM K-FILES(.02 TAPER)

size	21MM	25MM
#8	NT502-101	NT502-201
#10	NT502-102	NT502-202
#15	NT502-103	NT502-203
#20	NT502-104	NT502-204
#25	NT502-105	NT502-205
#30	NT502-106	NT502-206
#35	NT502-107	NT502-207
#40	NT502-108	NT502-208
#45	NT502-109	NT502-209
#50	NT502-110	NT502-210
#55	NT502-111	NT502-211
#60	NT502-112	NT502-212
#70	NT502-113	NT502-213
#80	NT502-114	NT502-214
#15/40	NT502-191	NT502-291
#45/80	NT502-192	NT502-292

(Pkg. of 6)

NITI H-FILES 
NICKEL-TITANIUM H-FILES(.02 TAPER)

size	21MM	25MM
#8	NT503-101	NT503-201
#10	NT503-102	NT503-202
#15	NT503-103	NT503-203
#20	NT503-104	NT503-204
#25	NT503-105	NT503-205
#30	NT503-106	NT503-206
#35	NT503-107	NT503-207
#40	NT503-108	NT503-208
#45	NT503-109	NT503-209
#50	NT503-110	NT503-210
#55	NT503-111	NT503-211
#60	NT503-112	NT503-212
#70	NT503-113	NT503-213
#80	NT503-114	NT503-214
#15/40	NT503-191	NT503-291
#45/80	NT503-192	NT503-292

(Pkg. of 6)

ENDODONTIC NICKEL TITANIUM HAND FILES/BLISTER PACKAGE



- Blister package for single use
- Increased Flexibility
- Strength
- Convenience
- Efficiency
- Safety

NITI REAMERS



NICKEL-TITANIUM REAMERS(.02 TAPER)

size	21MM	25MM
#8	NTB501-101	NTB501-201
#10	NTB501-102	NTB501-202
#15	NTB501-103	NTB501-203
#20	NTB501-104	NTB501-204
#25	NTB501-105	NTB501-205
#30	NTB501-106	NTB501-206
#35	NTB501-107	NTB501-207
#40	NTB501-108	NTB501-208
#45	NTB501-109	NTB501-209
#50	NTB501-110	NTB501-210
#55	NTB501-111	NTB501-211
#60	NTB501-112	NTB501-212
#70	NTB501-113	NTB501-213
#80	NTB501-114	NTB501-214
#15/40	NTB501-191	NTB501-291
#45/80	NTB501-192	NTB501-292

(Pkg. of 6)

NITI K-FILES



NICKEL-TITANIUM K-FILES(.02 TAPER)

size	21MM	25MM
#8	NTB502-101	NTB502-201
#10	NTB502-102	NTB502-202
#15	NTB502-103	NTB502-203
#20	NTB502-104	NTB502-204
#25	NTB502-105	NTB502-205
#30	NTB502-106	NTB502-206
#35	NTB502-107	NTB502-207
#40	NTB502-108	NTB502-208
#45	NTB502-109	NTB502-209
#50	NTB502-110	NTB502-210
#55	NTB502-111	NTB502-211
#60	NTB502-112	NTB502-212
#70	NTB502-113	NTB502-213
#80	NTB502-114	NTB502-214
#15/40	NTB502-191	NTB502-291
#45/80	NTB502-192	NTB502-292

(Pkg. of 6)

NITI H-FILES



NICKEL-TITANIUM H-FILES(.02 TAPER)

size	21MM	25MM
#8	NTB503-101	NTB503-201
#10	NTB503-102	NTB503-202
#15	NTB503-103	NTB503-203
#20	NTB503-104	NTB503-204
#25	NTB503-105	NTB503-205
#30	NTB503-106	NTB503-206
#35	NTB503-107	NTB503-207
#40	NTB503-108	NTB503-208
#45	NTB503-109	NTB503-209
#50	NTB503-110	NTB503-210
#55	NTB503-111	NTB503-211
#60	NTB503-112	NTB503-212
#70	NTB503-113	NTB503-213
#80	NTB503-114	NTB503-214
#15/40	NTB503-191	NTB503-291
#45/80	NTB503-192	NTB503-292

(Pkg. of 6)

ENDO STOP DISPENSER

• **ITEM #701-000**

- Contains 120 Assorted Round stops.
- Easy to use.
- Features ruler for measuring.



- Available Round Type.
- 100 Stops per poly bag.

White	701-101
Yellow	701-102
Red	701-103
Blue	701-104
Green	701-105
Black	701-106
Assorted	701-191