


DIA-ENDORING

• **ITEM #850-001**

- A convenient ring for Hand & Rotary Files.
- Gauges from 1 to 30mm.
- A simple & practical design for measuring files.
- Gold plated ring prevents color change.
- Dimension: 45 x 25 x 15m/m (W x H x D)



ENDO FINGERING

• **ITEM #790-101**

- Provides complete control of entire endodontic files.
- Disposable sponges hold files in a convenient position.
- Autoclavable except cups
- Contents :


- **ITEM #790-111 :** Disposable Cups(48pcs)
- **ITEM #790-121 :** Disposable Sponges(48pcs)
- Finger Rings(3pcs)
- Cup Rings(2pcs)



STERI ENDO GUARD W/RULER(ENGINE FILE)

• **ITEM #716-203 Yellow**
• **ITEM #716-206 Blue**


- Holds 16 NiTi Rotary Files.
- Right or left hand use.
- Available to insert sponge to clean or hold files.
- Diameter: 75m/m, Height: 35m/m



FILE MEASURE SET(PLASTIC)

• **ITEM #780-001**

- Holds 32 Measuring Grades.
- Gauge from 12 to 27.5m/m, Including half measurements.
- Dimensions: 87 x 31 x 20m/m (W x H x D)



DIAFLEX™ IRRIGATION SYRINGE FOR ROOT CANAL



• **ITEM #1003-301**

- Flexible soft polypropylene tip syringe for root canal irrigation.
- Autoclavable and chemically cleanable for multiple use up to 50 times.
- 1.3ml volume imprinted at intervals of 1mm.
- Provides maximum clinical effect in a minimum of time.
- DiaFlex is compatible with **Lure Lock Syringes**.


DiaFlex - III

- **ITEM #1003-301**
- Package of 3 color coded syringes, 20 changeable tips and silicone oil

DiaFlex - V

- **ITEM #1003-501**
- Package of 5 color coded syringes, 20 changeable tips and silicone oil

- **ITEM #1001-315**
- Blue 20 disposable tips in a **blister package**



Tip Specification