



STERI. IRRIGATION TIPS-SIDE OPEN TYPE

- Disposable Endodontic & Periodontal Irrigation Tips with Safety.
- Unique design with a side-port for superior effectiveness and patient comfort in irrigating root canals and periodontal pockets.
- Pre-sterilized by E.O Gas.
- Use with Lure Lock Syringes.
- Side-Port and Closed-end Irrigation tips.



PERIODONTAL POCKETS

- Side-port design make an effective irrigate pockets and closed, rounded tip can reduce the risk of tissue damage.



ROOT CANALS

- Side-port dispersal can create an upward flushing motion for the complete root canal preparations and cleanse the gingival sulci fully without damaging tissue.



size	Color	ITEM No.
31G (0.27 x 25mm)	Clear White	951-101
30G (0.30 x 25mm)	Yellow	951-102
28G (0.36 x 25mm)	Blue Green	951-103

(Pkg. of 50)



IRRIGATION TIPS-END OPEN TYPE



- Disposable Endodontic & Periodontal Irrigation tips
- Use with Lure Lock Syringes.
- Thoroughly irrigate root canal preparations.
- Round end tip can prevent the secondary injury during the use.
- Widest selection of tips: 6 sizes from 18G to 27G.

size	Color	ITEM No.
18G (1.2 x 25mm)	Pink	951-201
19G (1.1 x 25mm)	Cream	951-202
20G (0.9 x 25mm)	Yellow	951-203
22G (0.7 x 25mm)	Black	951-204
25G (0.5 x 25mm)	Orange	951-205
27G (0.4 x 25mm)	Medium Grey	951-206

(Pkg. of 100)

