


SCALING BRUSH LATCH TYPE

- **ITEM #941-001**
- Packed 72 brushes.
- Regular cup style.
- For polishing of tooth surface after scaling.
- Useful for wide space.
- Diameter: 4.49mm



SCALING POLISHING BRUSH LATCH TYPE

- **ITEM #941-103**
- Packed 36 brushes.
- Small cup style.
- Useful for narrow space.
- Diameter : 2.27mm



CANDLE BRUSH

- **ITEM #941-104 :**
Latch Type, Packed 24 brushes
- **ITEM #941-105 :**
Screw Type, Packed 24 brushes and 10 mandrels
- For removing stains from occlusal pits and fissures and for preparing tooth surface prior to sealant and orthodontic application.
- Diameter : Latch Type : 4.50mm / Screw Type : 5.00mm



POST BRUSH

- **ITEM #941-100 :**
Mandrels for Screw Type, packed 20 pcs
- **ITEM #941-101 :**
Latch Type, Packed 30 brushes
- **ITEM #941-102 :**
Screw Type, Packed 30 brushes
- For quick cleaning root canal preparation, under bridge.
- Available to insert sterilizing liquid during turn it left in Contra Angle.
- Available to flow out the stains like cement, liquid during turn it right.
- Recommended rotation: 300-800rpm.
- Single use.
- Diameter : Latch Type : 1.20mm / Screw Type : 1.03mm



PROPHY BRUSH

- **ITEM #950-101 :** Latch Type
- **ITEM #950-102 :** Screw Type
- Packed 144 brushes
- Diameter
Latch Type : 4.48mm
Screw Type : 4.48mm



PROPHY BRUSH POINTED HEAD

- **ITEM #950-111 :** Latch Type
- **ITEM #950-112 :** Screw Type
- Packed 144 brushes
- Diameter
Latch Type : 4.48mm
Screw Type : 4.48mm



PROPHY BRUSH LARGE HEAD LATCH TYPE

- **ITEM #950-121**
- Packed 144 brushes
- Diameter : 5.98mm



PROPHY BRUSH LARGE POINTED HEAD LATCH TYPE

- **ITEM #950-131**
- Packed 144 brushes
- Diameter : 5.98mm



SILICONE CUP FOR POLISHING LATCH TYPE

- **ITEM #941-002**
- Soft, white cups
- For polishing the tooth surface, achieves a high shine quickly.
- Use conveniently after insert in Contra Angle
- Pkg. of 72



PROPHY RUBBER CUP RIBBED & WEBBED

- **ITEM #950-201** : Latch Type
- **ITEM #950-202** : Screw Type
- Packed 144 cups



Latch Type (RC - 001)



Screw Type (RC - 002)

PROPHY RUBBER CUP WEBBED

- **ITEM #950-211** : Latch Type
- **ITEM #950-212** : Screw Type
- Packed 144 cups



Latch Type (RC - 003)



Screw Type (RC - 004)

PROPHY RUBBER CUP RIBBED

- **ITEM #950-221** : Latch Type
- **ITEM #950-222** : Screw Type
- Packed 144 cups



Latch Type (RC - 005)



Screw Type (RC - 006)

PROPHY RUBBER CUP UMBRELLA SHAPE

- **ITEM #950-231** : Latch Type
- **ITEM #950-232** : Screw Type
- Packed 144 cups



Latch Type (RC - 007)



Screw Type (RC - 008)

PROPHY RUBBER CUP SNAP-ON

- **ITEM #950-251** : Ribbed
- **ITEM #950-252** : Webbed
- Packed 144 cups



Ribbed (RC - 010)



Webbed (RC - 011)

PROPHY RUBBER CUP WEBBED WITH WASHER LATCH TYPE

- **ITEM #950-241**
- Packed 144 cups



(RC - 009)

PROPHY RUBBER CUP SNAP-ON RIBBED & WEBBED

- **ITEM #950-253**
- Packed 144 cups



(RC - 012)