

BUR BLOCK W/COVER



Type A

- ITEM # 836-001
-  
- Holds 24 Friction Grip(high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type B

- ITEM # 836-002
-  
- Holds 24 Right Angle(low speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type C

- ITEM # 836-003
-  
- Holds 12 Right Angle(low speed) Burs & 12 Friction Grip(high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type D

- ITEM # 836-004
-  
- Holds 36 Friction Grip(high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 50m/m (W x H x D)



Type E


- ITEM # 836-005
-  
- Holds 36 Right Angle(low speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 50m/m (W x H x D)

STERI BUR BLOCK W/COVER



Type A

- ITEM #837-001
-  
- Holds 24 Friction Grip(high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m(W x H x D)





Type B

- ITEM #837-002
-  
- Holds 24 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m(W x H x D)



Type C

- ITEM #837-003
-  
- Holds 12 Friction Grip(high speed) Burs & 12 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m(W x H x D)



Type D

- ITEM #837-004
-  
- Holds 36 Friction Grip(high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 50m/m(W x H x D)



Type E

- ITEM #837-005
-  
- Holds 36 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 50m/m(W x H x D)

MINI BUR STAND

Type A



- Silver ITEM #838-101
- Red ITEM #838-103
- Blue ITEM #838-104



Type B



- Silver ITEM #838-201
- Red ITEM #838-203
- Blue ITEM #838-204



Type C



- Silver ITEM #838-301
- Red ITEM #838-303
- Blue ITEM #838-304



- Holds 12 Friction Grip(high speed) Burs
- Durable Aluminum
- Dimensions: 45 x 32 x 25m/m (W x H x D)

- Holds 12 Right Angle(low speed) Burs
- Durable Aluminum
- Dimensions: 45 x 32 x 25m/m (W x H x D)

- Holds 9 Friction Grip(high speed) Burs & 3 Right Angle(low speed) Burs
- Durable Aluminum
- Dimensions: 45 x 32 x 25 m/m (W x H x D)



BUR BLOCK

Type A

- ITEM # 832-001



- Holds 48 Friction Grip (high speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)



Type B

- ITEM # 832-002



- Holds 48 Right Angle(low speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)



Type C

- ITEM # 832-003




- Holds 24 Friction Grip (high speed) Burs & 24 Right Angle(low speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)




STERI BUR GUARD(22 HOLE)

- ITEM #716-103 Yellow
- ITEM #716-106 Blue
- Fits 22 latch & friction grip burs
- Lid does not allow the burs to fall out of the stand.
- Diameter: 75m/m, Height: 35m/m



MAGNETIC BUR BLOCK(OBLONG-14 HOLE)

- Type A**
- ITEM #717-103 Yellow
 - ITEM #717-105 Beige
 - ITEM #717-106 Blue
 - Dimensions: 77 x 13 x 33m/m (W x H x D)



Type A

- Type B**
- ITEM #717-203 Yellow
 - ITEM #717-205 Beige
 - ITEM #717-206 Blue
 - Dimensions: 76 x 12 x 37mm (W x H x D)
 - 14 Holes
 - Fits both FG and RA types.



Type B

GERMICIDE ROUND BUR BOX

ITEM #915-001




- Contains perforated insert for drainage.
- Practical & Durable.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter: 75m/m, Height: 35m/m



AL BOX FOR BUR HOLDER(W/7 INSERTS)

- ITEM # 831-001
- Holds 24 High Speed Burs & 32 Low Speed Burs.
 - Interchangeable Inserts. Versatile.
 - Easy to transport.
 - Dimensions : 125 x 53 x 130m/m (W x H x D)



MAGNETIC BUR BLOCK(SQUARE-8 HOLE)

- ITEM #717-303 Yellow
- ITEM #717-306 Blue
- Fits both FG and RA types.
- Dimensions : 40 x 13 x 40m/m (W x H x D)



MAGNETIC BUR BLOCK(ROUND 24 HOLE)



- 8 Holes**
- ITEM #717-503 Yellow
 - ITEM #717-505 Beige
 - ITEM #717-506 Blue
 - Diameter : 43m/m, Height: 10m/m



8 Holes

- 24 Holes**
- ITEM #717-603 Yellow
 - ITEM #717-605 Beige
 - ITEM #717-606 Blue
 - Diameter : 73m/m, Height: 10m/m
 - Fits both FG & RA types



24 Holes

BUR STERI BATH

- ITEM # 730-001 White
- ITEM # 730-002 Transparent



- Used for Sterilizing Burs & Small Instruments.
- Removable insert for liquid drainage.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter : 88 m/m
- Height : 68 m/m



DOWEL PIN

- ITEM #907-001: Twin pin
- ITEM #907-002: Single pin



Twin pin with metal sleeve & Rubber cap

- Rubber cap saves working time.
- Rubber cap facilitates work.
- Exact fit into preparations and no distortion.
- Pkg.of 1,000

