

**GLASS DAPPEN DISHES**



- ITEM #601-001 Clear
- ITEM #601-002 Amber
- ITEM #601-003 Blue
- ITEM #601-004 Green



- Standard Dappen Dish.
- Do not Autoclave.
- Easy to clean.
- Sold Individually or Box of 12pcs.

**MEDICAMENT BOTTLES**



- Locking lids for good sealing capabilities.
- Available in 4 colors for coding of chemicals.
- Easy to use and store.
- High quality glass, economically priced.
- 10ml volume
- Do not Autoclave.

**ROUND TYPE**

- ITEM #604-001 Clear
- ITEM #604-002 Amber
- ITEM #604-003 Blue
- ITEM #604-004 Green

**SQUARE TYPE**

- ITEM #605-001 Clear
- ITEM #605-002 Amber
- ITEM #605-003 Blue
- ITEM #605-004 Green



**MIXING SLAB**



- High quality, competitively priced Mixing Slabs available in 5 sizes.
- Do not Autoclave.



- ITEM #603-005 LT  
- Dimensions : 150 x 15 x 80m/m (W x H x D)
- ITEM #603-004 LM  
- Dimensions : 150 x 10 x 75m/m (W x H x D)
- ITEM #603-003 M  
- Dimensions : 122 x 5 x 75m/m (W x H x D)
- ITEM #603-002 L  
- Dimensions : 112 x 6 x 76m/m (W x H x D)
- ITEM #603-001 S  
- Dimensions : 76 x 6 x 65m/m (W x H x D)

**CHAIR SIDE TABLE**



- ITEM #712-001 White
- Comes with attachments.
- Can be adapted to several types of chairs.
- Hole Diameter : available from 25 to 50m/m.
- Easy access to instruments.
- Dimensions : 390 x 30 x 235m/m (W x H x D)



**GLASS ENDO ORGANIZER**

- ITEM #602-000



- Shallow compartment walls for liquid sterilization.
- 6 inch diameter dish with 12 shallow wall compartments.
- Glass lid with depressions for chemicals.
- Do not Autoclave.



**MEDICAMENT BOTTLE STAND**

- ITEM #606-001



- Fits 5 Round or Square Medicament Bottles.
- Opaque high quality glass.
- Do not Autoclave.
- Color : Milky White
- Dimensions : 215 x 27 x 55m/m (W x H x D)



**BIOHAZARD HANDY HOOPER**

- ITEM #704-101 White



- A small plastic waste receptacle for biohazard items
- Can be wall mounted
- Dimensions : 170 x 160 x 100m/m(W x H x D)



**TISSUE HOLDER**

- ITEM #712-101 White



- Can be adapted to several types of chairs.
- Tissue box holder fits on unit chair pole.
- Hole Diameter : available from 25 to 50m/m.
- Dimensions : 245 x 54 x 125m/m(W x H x D)

